



Release Notes

MATRIXSPECS/: Validation (VALID)

Date: 2015/01/22 16:11:14

Document: MRM-1003-RN (0006)

Page: 1 / 2

RELEASE NOTES

SOFTWARE IDENTIFICATION

Matrix Requirements Medical 1.3.641.4838

NEW FUNCTIONALITY IN THIS RELEASE

Controlled Documents and Electronic Signatures

- Support for FDA CFR part 11

Time Zone definition and date / time formatting

- Setting of global time zone per and formatting server used as default and for all reports
- Overwrite for user interface per user

Password strength and ageing

- Possibility to require a certain password strength and
- Possibility to disable project access if password is too old

Possibility to upload project logos

Possibility to define a word template per project

Administration Client to manage

- Project access per user
- User identification, passwords, signatures etc
- Time zone and formatting default for server and per user overwrite
- Project logos
- JIRA setup
- Controlled document setup



Release Notes

MATRIXSPECS/: Validation (VALID)

Date: 2015/01/22 16:11:14

Document: MRM-1003-RN (0006)



Page: 2 / 2

KNOWN ISSUES

MATRIX-480 SIGN OUT and SIGN IN directly after does not work (Safari)

If a user does a sign out in safari, in order to sign in with a different user name, the sign in dialog will appear twice, since the safari 8 browser does not replace the invalid session cookie automatically.

SIGNATURES

Signature Meaning	Name	Title	Date	Signature
Approved By	Wolfgang Huber	CEO	2015/01/22	 sign-p14-i2307-b20150122161151553
Approved By	Yves Berquin	CEO	2015/01/22	 sign-p14-i2307-b20150122163648625